

## Semper 6 Command Reference

### qget

*This command is specific to...  
PC + Quantimet 520 greystore*

<b>keys:</b>	<b>[to]</b>	<i>&lt;number&gt;</i>	output Semper picture
	<b>name</b>	<i>'&lt;text&gt;'</i>	source image file (Quantimet 520 format)

Use **qget** to read a Quantimet 520 format image file into a Semper picture.

#### Examples

```
qget name '520PIC' to 2
```

This command reads the Quantimet 520 format disc file called *520PIC.IMG* into Semper picture 2.

```
qget 2 name '520PIC.IMG'
```

This command is identical in effect to the example given above.

#### Description

If you do not specify a filename using the **name** key, Semper prompts for a name. Note that the file extension defaults to *.IMG* but you can specify other extensions explicitly.

**qget** searches for files in the current directory and then throughout the DOS PATH. You can avoid a time-consuming path scan by specifying a full path name with the filename.

#### Notes

see also: **qput**

#### Defaults and Ranges

keys/options	defaults	range
<b>[to]</b>	current picture, held in the variable <i>select</i>	valid picture number
<b>name</b>	<i>none</i> ; if you do not specify a name, Semper prompts for a filename	valid filename